

Technical challenges

Jonas T. Larsson
SSI

Public Health WGS workshop


With focus on Illumina

- ✧ DNA preparation from cells
- ✧ Nextera XT Library prep
 - Short insert size
- ✧ Using different platforms, Illumina-Iontorrent-454-whatnot
 - HiSeq vs MiSeq
 - SNP pipeline
 - Gene by gene
- ✧ <http://seqanswers.com/>



How important is the DNA preparation?

- ❖ We have tried
 - Promega Wizard SV
 - Qiagen Gentra Puregene Yeast/Bact. Kit
- ❖ Our limited observation is that Nextera XT preps are identical from the same DNA prep
- ❖ Need to look closer into why some DNA preps seemingly have 'bad library properties'



DNA prep

Tagmentation

Reduced-cycle PCR

PCR purification


Quantify Libraries

Pool and run on
instrument

Short insert size



• Spotted in Illumina Sequencing Analysis Viewer



Short insert size

- ❖ Insert sizes with the default Nextera XT protocol are much too small
- ❖ Changing from 25 to 20 μ l AMPure XP beads in the PCR purification

DNA prep


Tagmentation

Reduced-cycle PCR

PCR purification

Quantify Libraries

Pool and run on instrument



Normalization

- ❖ We are now using a protocol from CDC without the bead-based normalization
- ❖ Single quantification with Qubit
- ❖ Simple dilution scheme
 - Compensate for genome sizes
- ❖ Skips microplate shaking/magnetic bead steps
- ❖ More even normalization

